

All Ceramic Materials for All-Ceramic Restorations

Zirtooth™


Zirconia Blocks for Dental CAD/CAM System



HASS BIO America, Inc.
HASS

Zirtooth™

Zirtooth™ Technology

Zirtooth™ zirconia blocks are made from 100% premium Tosoh  Yttria-stabilized tetragonal zirconia (YTZ) powders applied the latest hydrolysis process and nano-technology, resulting in higher flexural strength and fracture toughness.




A Shade

B Shade




Product Line-Up

Systems	Model No.	Dimensions (mm)	Zirtooth™ Multi	Zirtooth™ Fulluster
	S-4019	40 × 19 × 15.5	○	○
	S-5519	55 × 19 × 15.5	○	○
	S-6525	65 × 25 × 22	○	○

* Sirona is not the trademark of HASS Corp.

Available Shades

Model Name	Shades
Zirtooth™ Multi	A Type
	B Type
Zirtooth™ Fulluster	A0 / A1 / A2

hassbioamerica.com  hassbio
 Contact : 323-939-1668 | hassbioamerica@hassbio.com
 Technical Support : 213-458-2033

HASS BIO America, Inc.
 10400 Eaton Place, Suite 220, Fairfax, Virginia 22030

HASS BIO America, Inc.
HASS